

Chloride TC

UFI: 15YD-KXM8-951X-THSD



Danger

H314 - Causes severe skin burns and eye damage

H317 - May cause an allergic skin reaction

H319 - Causes serious eye irritation

H332 - Harmful if inhaled

H340 - May cause genetic defects

H350i - May cause cancer by inhalation

H360FD - May damage fertility. May damage the unborn child

H411 - Toxic to aquatic life with long lasting effects

EUH208 - Contains Potassium dichromate. May produce an allergic reaction

P201 - Obtain special instruction before use

P260 - Do not breathe dust/fume/gas/mist/vapours/spray

P273 - Avoid release to the environment

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P303 + P361 + P353 - IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing

P310 - Immediately call a POISON CENTER or doctor

P391 - Collect spillage

Test Kit

Cl⁻ (0 - 5000 mg/l)
NaCl (0 - 8000 mg/l)



Chloride TC



Chloride Tablet-count method



Chlorid Tabletten-Zählverfahren



Chlorure Methode de comptage à pilule



Cloruro Metodo di conteggio de pastiglie



Cloruro Método Contador de tabletas




Klorid Tabletællings-metoden



Chloride Tablet-telmethode

NOTRUF/EMERGENCY: +49 (0)89-19240 (24h)

Water-I.d. GmbH
Daimlerstr. 20
D-76344 Eggenstein
Tel. +49 (0) 721 - 782029-0
Fax +49 (0) 721 - 782029-11
www.water-id.com
info@water-id.com




0 - 100 mg/l Cl^-
(0 - 160 mg/l NaCl)

0 - 200 mg/l Cl^-
(0 - 320 mg/l NaCl)


0 - 1000 mg/l Cl^-
(0 - 1600 mg/l NaCl)

0 - 5000 mg/l Cl^-
(0 - 8000 mg/l NaCl)


1




2





3



4



5



NaCl

100 ml: (-1) x 16 = mg/l NaCl

50 ml: (-1) x 32 = mg/l NaCl

10 ml: (-1) x 160 = mg/l NaCl

2 ml: (-1) x 800 = mg/l NaCl

Cl^-

100 ml: (-1) x 10 = mg/l Cl^-

50 ml: (-1) x 20 = mg/l Cl^-

10 ml: (-1) x 100 = mg/l Cl^-

2 ml: (-1) x 500 = mg/l Cl^-

